



#### SPECIFICATIONS

Input voltage,	
Frequency range,	Input AC 100~240V~, 50/60Hz, 0.3A
Current	
Dimension	95mm x 125mm x 74mm (Body)
Weight	270g (Body)
Color	Rose Gold , Purple Black
Mode	GYM, SLIM
Application Area	Abdomen, Back, Forearms, Thighs, Calves



55, Teheran-ro 25-gil, Gangnam-gu, Seoul, Korea, 06126  
T. +82 2 3462 5400 F. +82 2 3462 5105 E. [gtg@gtgwellness.co.kr](mailto:gtg@gtgwellness.co.kr)

[www.gtgwellness.co.kr](http://www.gtgwellness.co.kr)

## LEBODY — FORM —



# Discover. Your Body Line.





## MID-FREQUENCY. FOR YOUR BODY.

LEBODY FORM helps to improve skin elasticity and swollen body by contraction and relaxation of muscles through MID-FREQUENCY.

Cellulite reduction

Decrease body waste and body fat

Improve skin elasticity and skin repair effect

Strengthen muscles

Part used



## PATENTED BY 'SEPPIC', FRANCE.

Adipoless & Herbal blended Botaniceutical plus-10.  
Doubled effectiveness for slim, elastic body.



## LEBODY FORM PATENT

Mid-frequency abdominal fat reduction device

KOREA-Patent No. 10-1591736

Portable mid-and low-frequency device capable of multi directional muscle movement

KOREA-Patent No. 10-1646086

PCT / KR2016 / 009164



## 2 MODES, HARDER WITH BEAT FREQUENCY.

Manage daily by Gym mode, better after work out. Intense body care by Slim mode, more muscle contraction by BEAT FREQUENCY.



Illuminance Sensor  
Electrode  
Aesthetic LED

## SAFE AND SECURE.

Safe use by 6 step intensity control and failsafe sensor.

## EXPERIENCE THE EFFECTS. CLINICALLY PROVEN.

Possesses most abdominal researches. Authorized and verified from official agencies, with high satisfaction.

- ✓ Temporarily Reduces Cellulite
- ✓ Weakened Abdominal Internal Elasticity Improvements
- ✓ Unbalanced Abdominal Skin Texture Improvements
- ✓ Abdominal and Waist Swelling Improvements
- ✓ Upper and Lower Abdominal Swelling Improvements
- ✓ Thigh Swelling Improvements
- ✓ Weakened Forearm Elasticity Improvements
- ✓ Weakened Thigh Elasticity Improvements
- ✓ Completing Skin Irritation Test

